|   |                         |    |
|---|-------------------------|----|
| Ⓔ | IPL Hair Removal Device | 1  |
| Ⓕ | IPL-Haarentferner       | 12 |
| Ⓖ | IPL Hårborttagare       | 22 |
| Ⓛ | IPL fotoepiliatorius    | 33 |
| Ⓛ | Fotoepilators IPL       | 44 |
| Ⓔ | Fotoepilaator IPL       | 56 |

### **Contents**

#### **Welcome**

#### **In the package you will find**

#### **How IPL works**

#### **Hair and skin tones**

#### **What to expect**

#### **Who should not use Hilda? Contraindication**

- General conditions

- Medications/History

- Pathologies/Disorders

- Skin condition

- Location/areas

#### **Important**

- Danger

- Warning

- To prevent damage

- Caution

#### **How to use**

- Skin test

- Selecting the right light intensity

- Recommended treatment schedule

- Handling the device

- Two treatment modes: Stamp & Flash and Slide & Flash

#### **Aftercare**

- Common skin reactions

- Rare side effects

- Aftercare

#### **Using your Hilda before and after tanning**

- Tanning with natural or artificial sunlight

- Tanning with creams

#### **Recycling**

#### **Technical specifications**

#### **Troubleshooting**

#### **Warranty and support**

# Hilda

## Welcome

Welcome to the beauty world of Hilda! Hilda uses Intense Pulsed Light (IPL) technology, known as one of the most effective methods to continuously prevent hair regrowth. Hilda is gentle and offers convenient and effective treatment at a light intensity that you find comfortable. Unwanted hairs are finally a thing of the past. Enjoy the feeling of being hair-free and look and feel amazing every day.

## In the package you will find

- IPL hair removal device;
- power cord;
- instructions;
- safety glasses.

## How IPL works

With Intense Pulsed Light technology, gentle pulses of light are applied to the skin and absorbed by the hair root. The lighter the skin and the darker the hair, the better the pulses of light are absorbed.

The pulses of light stimulate the hair follicle to go into a resting phase. As a consequence, the hair sheds naturally and hair regrowth is prevented.

## Hair and skin tones

The intense pulsed light device works best for removing black, dark brown, medium or light brown hair. Intense pulsed light technology cannot affect light, red or grey hair. Also, hair removal procedures are not effective for people with very dark skin.

### Suitable skin tones

### Suitable hair shades

## What to expect

- After the first treatment, it can take 2 to 4 weeks for the hairs to fall out. In the first weeks following the initial treatments, you still see some hairs growing. These are likely to be hairs that were not in their growing phase during the first treatments.
- After 4-6 treatments, you should see a noticeable reduction in hair growth. However, to effectively treat all hairs, it is important to keep on treating according to the recommended treatment schedule.
- After 5-7 treatments, you should see a significant reduction of hair growth in the areas that you treated with Hilda. A reduction of hair density should be visible as well.

## Who should not use Hilda? Contraindication

### General conditions

- Never use the device if you have skin type VI (You rarely ever sunburn, very dark tanning). In this case you run a high risk of developing skin reactions, such as hyperpigmentation and hypopigmentation, strong redness or burns.
- Never use the device if you are pregnant or breastfeeding as the device was not tested on pregnant or breastfeeding women.
- Never use the device if you have any active implants such as a pacemaker, neurostimulator, insulin pump etc.

### Medications/History

Never use the device if you take any of the medications listed below:

- If your skin is currently being treated with or has recently been treated in the past week with Alpha-Hydroxy Acids (AHAs), Beta-Hydroxy Acids (BHAs), topical isotretinoin and azelaic acid.
- If you have taken any form of isotretinoin Accutane or Roaccutane in the last six months. This treatment can make skin more susceptible to tears, wounds and irritations.
- If you are taking photosensitizing agents or medications, check the package insert of your medicine and never use the device if it is stated that it can cause photo-allergic reactions, photo-toxic reactions or if you have to avoid sun when taking this medicine.
- If you take anticoagulation medications, including heavy use of aspirin, in a manner which does not allow for a minimum 1-week washout period prior to each treatment.

Never use the device:

- If you have received radiation therapy or chemotherapy within the past 3 months.
- If you are on painkillers which reduce the skin's sensitivity to heat.
- If you take immunosuppressive medications.
- If you have had surgery in the areas to be treated in the last 3 weeks.

### Pathologies/Disorders

Never use the device:

- If you have diabetes or other systemic or metabolic diseases.
- If you have congestive heart disease.
- If you have a disease related to photosensitivity, such as polymorphic light eruption (PMLE), solar urticaria, porphyria etc.
- If you have a history of collagen disorder, including a history of keloid scar formation or a history of poor wound healing.
- If you have epilepsy with flashlight sensitivity.
- If your skin is sensitive to light and easily develops a rash or an allergic reaction.
- If you have a skin disease such as active skin cancer, you have a history of skin cancer or any other localised cancer in the areas to be treated.
- If you have a history of vascular disorder, such as the presence of varicose veins or vascular ectasia in the areas to be treated.
- If you have any bleeding disorders.
- If you have a history of immunosuppressive disease (including HIV infection or AIDS).

### Skin condition

Never use the device:

- If you have infections, eczema, burns, inflammation of hair follicles, open lacerations, abrasions, herpes simplex (cold sores), wounds or lesions and haematomas in the areas to be treated.
- On irritated (red or cut), sunburned, recently tanned or fake- tanned skin.

- On the following areas without consulting your doctor first: moles, freckles, large veins, darker pigmented areas, scars and skin anomalies. This can result in a burn and a change in skin colour, which makes it potentially harder to identify skin-related diseases.
- On following areas: warts, tattoos or permanent make-up.

## **Location/areas**

Never use the device on the following areas:

- Around the eyes and on or near the eyebrows.
- On lips, nipples, areolas, labia minora, vagina, anus and the inside of the nostrils and ears.
- Men must not use the device on the face and neck including all beard-growing areas, nor on the whole genital area.
- In areas where you use long-lasting deodorants. This can result in skin reactions.
- Over or near anything artificial like silicone implants, subcutaneous injection ports (for instance an insulin dispenser) or piercings.

Note: This list is not exhaustive. If you are not sure whether you can use the device, we advise you to consult your doctor.

## **Important**

### **Danger**

- Keep the device and the adapter dry.
- If the device is broken, do not touch any inner part to avoid electric shock.
- Water and electricity are a dangerous combination. Do not use this appliance in wet surroundings (e.g. near a filled bath, a running shower or a filled swimming pool).

### **Warning**

- This device is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, unless they have been given supervision or instruction concerning use of the device by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the device.
- The device is not intended for children under the age of 15 years. Teenagers aged between 15 and 18 years can use the device with the consent and/or assistance of their parents or the persons who have parental authority over them. Adults of 18 years and older can use the device freely.
- Always check the device before you use it. Do not use the device or adapter if it is damaged. Always replace a damaged part with one of the original type.
- Do not use the device if the UV filter of the light exit window is broken.
- Do not modify or cut off any part of the adapter or the cord, as this causes a hazardous situation.
- Do not use any pencil or pen to mark the areas to be treated. This may cause burns on your skin.
- If you have dark-coloured skin, be careful with treating a darker area immediately after a lighter area. The skin tone sensor may not immediately block the treatment on the darker body area.
- Hair removal by intense pulsed light sources can cause increased hair growth in some individuals. Based upon currently available data, the highest risk groups for this response are females of Mediterranean, Middle Eastern and South Asian heritage treated on the face and neck.
- Adapter, light exit window and the filter of the attachments can become very hot after usage. Do not touch the adapter, inner part of the light exit window and the filter or the inner part of the attachments without having these cool down.
- This device has a detachable power adapter. Only use the adapter provided with this device.

Note: If you notice a skin tone change from last treatment (e.g. due to tanning), we recommend you to perform a skin test and to wait 30 minutes before your next treatment.

## To prevent damage

- Make sure that nothing obstructs the airflow through the ventilation slots of the device.
- Never subject the device to heavy shocks and do not shake or drop it.
- If you take the device from a very cold environment to a very warm environment or vice versa, wait approximately 3 hours before you use it.
- Store the device in a dust free and dry place.
- Do not expose the device to temperatures lower than 15°C or higher than 35°C during use.
- To prevent damage, do not expose the device to direct sunlight or UV light for several hours.

## Caution

- This device is only intended for removing unwanted body hair from areas below the cheekbones. Do not use it for any other purpose. Doing so may expose you to a hazardous situation. Men must not use it on the face and neck including all beard-growing areas and the whole genital area.
- This device is not washable. Never immerse the device in water and do not rinse it under the tap.
- For hygienic reasons, the device should only be used by one person.
- Use the device only at settings suitable for your skin type. Use at higher settings than those recommended can increase the risk of skin reactions and side effects.
- Never use compressed air, scouring pads, abrasive cleaning agents or aggressive liquids such as petrol or acetone to clean the device.
- Do not look at the flash while using the device. Wear safety glasses and do not remove them while using the device. Use the device in a well-lit room so that the light is less glaring to your eyes.
- Never leave the device unattended when it is switched on. Always switch off the device after use.
- Do not use the device if any of the conditions mentioned in chapter 'Who should not use Hilda? Contraindication applies to you.
- Tanning with natural or artificial sunlight might influence the sensitivity and colour of your skin. Perform a skin test to determine the appropriate light intensity setting.
- Before you use Hilda, you should clean your skin and make sure it is hair-free, entirely dry and free from oily substances.
- Do not treat the same skin area more than once during a session. This does not improve the effectiveness of the treatment, but increases the risk of skin reactions.
- Hilda should never be painful. If you experience discomfort, reduce the light intensity setting.

## How to use

### Skin test

When you use Hilda for the first time or after recent tanning, perform a skin test on each area to be treated. The skin test is necessary to check your skin's reaction to the treatment and to determine the correct light intensity setting for each body area.

1. Do not try the device on difficult or sensitive areas (ankle and bony area).
2. Turn on the device.
3. Place the device at a 90° angle on the skin so that the integrated safety system is in contact with your skin. The integrated safety system prevents unintentional flashing without skin contact.
4. Press the flash button to release a flash.
5. Slide the device over the skin to the next area to be treated.
6. Increase the setting by one level, apply a flash and slide the device to the next area. Repeat this for all levels, as long as each level still feels comfortable.
7. After the skin test, wait 24 hours and check your skin for any reaction. If your skin shows reactions, choose the highest setting that did not result in any skin reaction for subsequent use.

## Selecting the right light intensity

Hilda provides 5 different light intensities and advises you on the right setting you should use based on your skin tone. You will always be able to change to a light intensity setting you find comfortable.

Skin types:

1. Very white, milk skin
2. Light skin
3. Beige skin
4. Brown skin
5. Dark brown skin
6. Black skin

| Skin types                                             | 1                     | 2                     | 3                     | 4                     | 5                     | 6 |
|--------------------------------------------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|---|
| Test (48 hours before hair removal for the first time) | 1 Intensity level     | 1 Intensity level     | 1 Intensity level     | 1 Intensity level     | 1 Intensity level     | x |
| First procedure                                        | 3 Intensity level     | 3 Intensity level     | 2 Intensity level     | 2 Intensity level     | 1 Intensity level     | x |
| Later hair removal                                     | 4 > 5 Intensity level | 4 > 5 Intensity level | 3 > 4 Intensity level | 3 > 4 Intensity level | 2 > 3 Intensity level | x |

Press the on/off button to turn on the device. The device starts to operate at light intensity setting 1.

To adjust the light intensity setting manually, use the toggle button until you have reached the required setting.

Using Hilda should never be painful. If you experience discomfort, reduce the light intensity setting. You can do this by using the up-down buttons.

## Recommended treatment schedule

The full cycle of procedures is 8 weeks.

Week 1-2: 2 treatments per week.

Week 3-4: 1 treatment per week.

Week 5-8: 1 treatment every 2 weeks.

Once the desired result is achieved, it is necessary to continue using the device once every 1-3 months for further inhibition of hair growth. Individual results may vary.

## Handling the device

1. Shave the hair and clean the skin dry.
2. Put the small plug in the device and put the adapter in the wall socket.
3. Switch on the device and select the right light intensity for your skin tone. For selecting the right skin tone, you can make use of the setting advice, see chapter: 'Selecting the right light intensity'.
4. Wear safety glasses and do not remove them while using the device.
5. Place the device at a 90° angle on the skin so that the integrated safety system is in contact with your skin. The integrated safety system prevents unintentional flashing without skin contact.
6. Press the device firmly onto your skin to ensure proper skin contact.
7. Press the flash button to release a flash. The device makes a soft popping sound. You should feel a warm sensation because of the flash. Do not look at the flash while using the device. Use the device in a well-lit room so that the light is less glaring to your eyes.

8. Place the device on the next area to be treated. After each flash, it takes up to 2 seconds until the device is ready to flash again.
9. To make sure that you have treated all areas, release the flashes close to each other. Effective light only comes out of the light exit window. Part of the device that is in contact with your skin is slightly larger therefore there should be some overlap. However, make sure you flash the same area only once.

Note: Flashing the same area twice does not improve the effectiveness of the treatment, but increases the risk of skin reactions.

Note: Do not use any pencil or pen to mark the areas to be treated, as this may cause side effects to your skin.

10. When you have finished the treatment, press and hold the on/off button to switch off the device. Remove the adapter from the wall socket.
11. After the hair removal procedure, it is recommended to apply a moisturiser.

## **Two treatment modes: Stamp & Flash and Slide & Flash**

Your Hilda has two treatment modes for more convenient use on different body areas:

- Stamp mode - performs one flash each time the button is pressed. It is designed for body parts that require gentle care and accuracy, such as the armpits, bikini area and face.
  - Automatic mode – aimed at removing hair from a large area, such as legs or arms. When this mode is activated, the device flashes automatically.
- With a short press of a button, you make a single flash. Press and hold to activate the automatic mode.

Note: The device needs up to 2 seconds in between two flashes, so make sure you flash on every area you wish to treat. Sliding too fast over the skin will cause missed spots.

## **After use**

### **Common skin reactions**

Your skin may show slight redness and/or may prickle, tingle or feel warm. This reaction is absolutely harmless and disappears quickly. Dry skin and itching may occur due to shaving or a combination of shaving and light treatment. You can cool the area with an ice pack or a wet face cloth. If dryness persists, you can apply a non-scented moisturiser on the treated area.

### **Rare side effects**

- Burns, excessive redness (e.g. around hair follicles) and swelling: these reactions occur rarely. They are the result of using a light intensity that is too high for your skin tone. If these reactions do not disappear within 3 days, we advise you to consult a doctor. Wait with the next treatment until the skin has healed completely and make sure you use a lower light intensity.
- Skin discoloration: this occurs very rarely. Skin discoloration manifests itself as either a darker patch (hyperpigmentation) or a lighter patch (hypopigmentation) than the surrounding area. This is the result of using a light intensity that is too high for your skin tone. If the discoloration does not disappear within 2 weeks, we advise you to consult a doctor. Do not treat discoloured areas until the discoloration has disappeared and your skin has regained its normal skin tone.
- Skin infection is very rare but is a possible risk following a (micro)wound, a skin burn, skin irritation etc.
- Epidermal heating (a sharply defined brownish area which often occurs with darker skin tones and is not accompanied with skin dryness): This reaction occurs very rarely. In case this reaction does not disappear within 1 week, we advise you to consult a doctor. Wait with the next treatment until the skin has healed completely and make sure you use a lower light intensity.
- Blistering (looks like small bubbles on the surface of the skin): this occurs very rarely. In case this reaction does not disappear within 1 month or when the skin gets infected, we advise you to consult a doctor. Wait with the next treatment until the skin has healed completely and make sure you use a lower light intensity.
- Scarring: often the secondary effect of a burn, which can take longer than a month to heal.

- Folliculitis (swelling around hair follicles combined with pustule formation): this reaction occurs very rarely and is the result of bacteria penetrating the damaged skin. In case this reaction occurs, we advise you to consult a doctor as folliculitis may need antibiotic ointment.

- Excessive pain: this can occur during or after treatment if you have used the device on skin that is not hair-free, if you use the device at a light intensity that is too high for your skin tone, if you flash the same area more than once and if you use the device on open wounds, inflammations, infections, tattoos, burns, etc.

## Aftercare

After use, you can safely apply lotions, creams, deodorant, moisturiser or cosmetics to the treated areas. If you experience skin irritation or skin redness after treatment, wait until it disappears before applying any product to your skin. If you experience skin irritation after applying a product to your skin, wash it off with water.

Note: when following all instructions you may still get skin reactions. In this case stop using the device and contact the consumer care centre in your country.

## Using your Hilda before and after tanning

### Tanning with natural or artificial sunlight

Intentionally exposing your skin to natural or artificial sunlight with the aim of developing a tan influences the sensitivity and colour of your skin. Therefore the following is important:

- After each treatment, wait at least 48 hours before tanning. Even after 48 hours, make sure that the treated skin does not show any redness from the treatment anymore.

- In case of exposing your skin to the sun (without tanning intentionally) in the 48 hours after treatment, use a sunblock SPF 50+ on the treated areas. After this period, you can use a sunblock SPF 30+ for two weeks.

- After tanning, wait at least 2 weeks before you use Hilda.

- After recent tanning, perform a skin test to determine the appropriate light intensity setting. For instructions see chapter 'Skin test'.

- Do not use Hilda on sunburned body areas.

Note: Occasional and indirect sun exposure does not qualify as tanning.

### Tanning with creams

If you have used an artificial tanning lotion, wait until the artificial tan has disappeared completely before you use the device.

## Recycling

- This symbol means that this product shall not be disposed of with normal household waste (2012/19/EU).

Follow your country's rules for the separate collection of electrical and electronic products and rechargeable batteries. Correct disposal helps prevent negative consequences for the environment and human health.

## Technical specifications

- Weight of the device: 0.3 kg;
- Total package weight: 1.23 kg;
- Device size: 140 \* 200 \* 40 mm;
- Package size: 260 \* 330 \* 67 mm;
- Voltage: 100v-240v, 12V;
- Power supply: AC
- Storage conditions temperature: -25° to +70°C;
- Storage conditions humidity: Less than 90% (non condensing);
- Emitted wavelengths: 620 nm

## Troubleshooting

| Problem                                                                                                   | Possible cause                                                                                                         | Solution                                                                                                                                                                                                                                                                                             |
|-----------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The device/adapter becomes warm during use.                                                               | It is normal for the device and adapter to become warm (but not too hot to touch) during use.                          | Use the device in a slightly cooler environment and/ or let it cool down before continuing use.                                                                                                                                                                                                      |
| The ready to flash light does not light up white.                                                         | The device is not completely in contact with your skin.                                                                | Place the device at a 90 degree angle on the skin so that the integrated safety system is in contact with your skin.                                                                                                                                                                                 |
| The device produces a strange smell.                                                                      | You have not removed the hairs on the area to be treated properly. These hairs may get burned and can cause the smell. | Pretreat your skin before you use Hilda.                                                                                                                                                                                                                                                             |
| The skin feels more sensitive than usual during treatment. I experience discomfort when I use the device. | The light intensity setting you use is too high.                                                                       | Check if you have selected the right light intensity setting. If necessary, select a lower setting.                                                                                                                                                                                                  |
|                                                                                                           | You did not remove the hairs on the areas to be treated.                                                               | Pretreat your skin before you use Hilda.                                                                                                                                                                                                                                                             |
|                                                                                                           | You treated an area for which the device is not intended.                                                              | Never use the device on the following areas: inner labia, anus, nipples, areolas, lips, moles, freckles, tattoos, piercings. inside the nostrils and ears, around the eyes and near the eyebrows. Men must not use it on the face and neck including all beard growing areas and whole genital area. |
| The skin reaction after treatment lasts longer than usual.                                                | You have used a light intensity setting which is too high for you.                                                     | Select a lower intensity next time.                                                                                                                                                                                                                                                                  |
| The results of the treatment are not satisfactory.                                                        | You have used a light intensity setting which is too low for you.                                                      | Select a higher setting next time.                                                                                                                                                                                                                                                                   |
|                                                                                                           | You did not flash an area adjacent to an area you treated before.                                                      | You have to release the flashes close to each other and there should be some overlap between the flashes.                                                                                                                                                                                            |
|                                                                                                           | The device is not effective on your body hair colour.                                                                  | If you have light blond, grey, red or white hair, the treatment is not effective.                                                                                                                                                                                                                    |
|                                                                                                           | You do not use the device as often as recommended.                                                                     | To remove all hairs successfully, we advise you to follow the recommended treatment schedule.                                                                                                                                                                                                        |

|  |                                           |                                                                                                                      |
|--|-------------------------------------------|----------------------------------------------------------------------------------------------------------------------|
|  |                                           | You can reduce the time between treatments, but do not treat more often than once every two weeks.                   |
|  | You respond more slowly to IPL treatment. | Continue using the device for at least 6 months, as hair regrowth can still decrease over the course of this period. |

**Manufacturer:**

Shenzhen youjuzhineng Technology Co., Ltd  
 503, No.2 Building, Songbai Industrial Zone, No. 4,  
 Yangyong Industrial Road, Tangxiayong Community,  
 Yanluo Street, Bao'an District, Shenzhen  
 Tel: +86 15986769873  
 Email: 736141031@qq.com

## Warranty

*This sheet contains information about the warranty conditions of the product. For detailed information and instructions on the use and maintenance of the device, see the instructions included with the order.*

### Warranty claim form

|                        |  |
|------------------------|--|
| Product name           |  |
| Date of purchase       |  |
| Business name, address |  |

#### 1. Warranty period

The warranty is valid from the date the purchaser purchased this product from the seller. The warranty period is 24 months.

#### 2. Warranty conditions

The warranty is valid upon presentation of the original product purchase document.

If the product is defective during the warranty period for any reason not covered by the warranty, the manufacturer undertakes to repair or replace the product. The warranty does not cover regular inspection, maintenance and repair or replacement of normally worn parts.

The warranty does not cover damage or malfunctions caused by:

- Improper use of the product and / or procedures not in accordance with the manufacturer's instructions for use, technical requirements or safety standards, including improper maintenance, falling to the ground and severe shocks.
- corrosion, dirt, water.
- Repairs performed outside of warranty service.
- Improper cleaning of the appliance, notwithstanding the recommendations given.
- Accidents or incidents or any causes beyond the control of the manufacturer.
- If the product does not meet the standards or specifications of a specific country other than the country of purchase.

Any attempt to adapt the product to the technical or safety requirements of other countries will void all warranty service. If no defects are found in the product, all service-related costs are paid by the buyer.

The warranty claim must be submitted by e-mail within two weeks of the defect occurring, describing in detail the basis for the warranty claim. In the event that the application is recognized as valid, the guarantor may, at its discretion:

- replace the defective product with a defective one,
- free of charge to find any defects found,
- to pay financial compensation when the elimination of defects will not be possible or if the elimination of the defect requires excessive costs depending on the value of the defect-free product and the type or significance of the defect found.

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